



INFORMATION FOR USE (IFU)

MDXR0027

Document #: MDX-SAL-IFU-017

Revision No: 001

Effective Date: 12-5-2025

Mock Stool Diluent

Description: Mock Stool Diluent is formulated to mimic a normal clinical stool sample.

Intended Use: Mock Stool Diluent may be used for research and development, analytical testing, and/or *in vitro* manufacturing of your diagnostic test.

Product Specifications:

| MRNDx Part Number | Volume | Storage Conditions | Stability |
|-------------------|--------|--------------------|-----------|
| MDXR0027-100 | 100mL | -20 to -80 °C | 24 months |
| MDXR0027-250 | 250mL | -20 to -80 °C | 24 months |
| MDXR0027-500 | 500mL | -20 to -80 °C | 24 months |
| MDXR0027-1000 | 1000mL | -20 to -80 °C | 24 months |

For Research Use Only – This product is intended for research and development use and should not be used for clinical diagnostic purposes.

IMPORTANT: Biosafety Level 1: The materials used to formulate this product contain no infectious material that pose hazard to laboratory personnel and the environment. These products are NOT intended for use in the manufacturing and/or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

This product was manufactured in a facility under a Quality Management System certified by ISO 13485:2016 for the manufacturing of medical devices.