



INFORMATION FOR USE (IFU)

MDXB0099

Document #: MDX-SAL-IFU-014

Revision No: 002

Effective Date: Captured Electronically

Normal Pooled Human Plasma

Description: Normal Human Plasma Pool is derived from whole blood collected from consented, healthy donors in an AABB–certified, IRB-approved, and FDA-licensed blood center. All units are screened negative for FDA required viral markers using FDA-cleared test methods. Various anticoagulants are available for collection, which include: K2 EDTA, K3 EDTA, CPD, Lithium Heparin, Sodium Heparin, Potassium Oxalate, and more.

Intended Use: Normal Human Plasma Pool may be used for research and development, analytical testing, and/or for *in vitro* manufacturing of your diagnostic test.

Product Specifications:

MRNDx Part Number	Volume	Storage Conditions	Stability
MDXB0099-10	10mL	-10 to -30 °C	24 months
MDXB0099-50	50mL	-10 to -30 °C	24 months
MDXB0099-100	100mL	-10 to -30 °C	24 months
MDXB0099-500	500mL	-10 to -30 °C	24 months

For Research Use Only – This product is intended for research and development use and should not be used for clinical diagnostic purposes.

Caution Human Biological Material: Biosafety Level 2: Includes moderate-risk agents that are present in the community and associated with human. All units are screened negative for FDA required viral markers using FDA-cleared test methods, however, all materials should be considered potentially infectious and handled with universal precautions.

IMPORTANT: These products are NOT intended for use in the manufacturing and/or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

This product was manufactured in a facility under a Quality Management System certified by ISO 13485:2016 for the manufacturing of medical devices.